

## Differences between bright-field, dark-field and phase contrast

Different imaging techniques are available for examining different specimens.

	BRIGHT-FIELD	DARK-FIELD	PHASE CONTRAST
<b>Principle of imaging</b>	<ul style="list-style-type: none"> <li>▪ Different levels of light absorption on different objects in the specimen</li> </ul>	<ul style="list-style-type: none"> <li>▪ Deflection of light on objects in the specimen</li> </ul>	<ul style="list-style-type: none"> <li>▪ Phase change when radiating through objects</li> </ul>
<b>Suitable for which samples</b>	<ul style="list-style-type: none"> <li>▪ High contrast samples</li> <li>▪ Coloured samples</li> <li>▪ Perfect for the "first look" at the sample</li> </ul>	<ul style="list-style-type: none"> <li>▪ Low contrast samples</li> <li>▪ Non-coloured samples</li> </ul>	<ul style="list-style-type: none"> <li>▪ Very thin biological specimens</li> <li>▪ Low contrast samples</li> <li>▪ Non-coloured samples</li> <li>▪ Living objects</li> </ul>
<b>Advantages</b>	<ul style="list-style-type: none"> <li>▪ Very easy</li> <li>▪ Very fast</li> <li>▪ Flat surface structures clearly visible</li> <li>▪ Correct colour impression</li> <li>▪ Components available in almost every optical microscope</li> </ul>	<ul style="list-style-type: none"> <li>▪ Simple</li> <li>▪ Samples with low contrast in bright field or almost transparent samples can be observed very well</li> <li>▪ Elevations on objects are easier to see than in bright-field</li> </ul>	<ul style="list-style-type: none"> <li>▪ Samples with low contrast in bright field or almost transparent samples can be observed very well with phase contrast</li> <li>▪ Elevations on objects are easier to see than in bright-field</li> <li>▪ Flat surface structures are clearly visible in contrast to dark-field</li> </ul>
<b>Disadvantages</b>	<ul style="list-style-type: none"> <li>▪ Low contrast in many samples, especially biological samples</li> <li>▪ Almost transparent samples are barely visible</li> <li>▪ Elevations on objects are difficult to recognise</li> </ul>	<ul style="list-style-type: none"> <li>▪ Not suitable for thick specimens</li> <li>▪ Flat surface structures of objects are poorly or not at all recognisable</li> <li>▪ Wrong colour impression</li> <li>▪ Special condenser required</li> <li>▪ High light intensity required, which can damage samples</li> <li>▪ Contamination is very clearly visible</li> </ul>	<ul style="list-style-type: none"> <li>▪ Not suitable for thick and medium-thick specimens</li> <li>▪ Complex adjustment of the phase contrast device on the microscope</li> <li>▪ Special phase contrast device and phase contrast objectives required</li> <li>▪ Phase contrast objectives lead to a loss of contrast, resolution and colour when used in bright-field</li> </ul>