

Troubleshooting MBL4000 series

| FAULT | POSSIBLE CAUSE | REMEDY |
|--|---|--|
| No microscope image visible | • Objective not correctly engaged in the revolving nosepiece | → Turn the revolving nosepiece until it clicks into place |
| | • Condenser not correctly inserted into the condenser mount | → Push the condenser into the condenser mount as far as possible |
| | • Illumination not switched on or defective | → Turn on the illumination → Contact customer service |
| | • Microscopy is started with very high magnification and specimen is not found | → Always start microscopy with the smallest magnification |
| | • A cover glass that is too thick was used, which means that the image cannot be focused without collision | → Use cover glasses with a thickness of 0.17 mm |
| Microscopy image is blurred | • Objective has collided with specimen | → Prevent objective from colliding with specimen |
| | • The 100X objective is used in air | → Always operate 100X objective with immersion oil |
| | • A cover glass that is too thick was used, which means that the image cannot be focused without collision | → Use cover glasses with a thickness of 0.17 mm |
| | • No cover glass was used | → Use cover glasses with a thickness of 0.17 mm |
| | • The specimen between the slide and the cover glass is thicker than necessary | → Microscopy with only as much specimen as necessary |
| | • Dioptre compensation is not set correctly, or dioptre compensation is incorrectly used for focussing the microscope image | → Do not use dioptre compensation to focus the microscope image → Perform dioptre compensation |
| There is dirt visible in the microscope image | • There is dirt on the slide or cover glass | → Use only clean slides and cover glasses and avoid fingerprints. |
| | • Dirt is on or in the eyepiece | → Carefully clean the eyepiece lens from both sides. To do this, use compressed air or a fine, lint-free cloth and water/Isopropanol. |
| | • There is dirt on or in the objective, e.g. immersion oil | → Carefully clean the lens of the objective. To do this, use compressed air or a fine, lint-free cloth and water/Isopropanol. → If there is dirt inside the objective, carefully(!) unscrew it and clean the lenses with extreme caution. |
| | • Dirt is in the remaining optics in the microscope head | → Remove the microscope head, remove the eyepieces and blow through the microscope head with compressed air |

| FAULT | POSSIBLE CAUSE | REMEDY |
|--|--|---|
| Contrast too weak | • Aperture diaphragm not set correctly | ➔ Set aperture diaphragm |
| | • The specimen between the slide and the cover glass is thicker than necessary | ➔ Microscopy with only as much specimen as necessary |
| | • The wrong aperture was set in the phase contrast device | ➔ Select the correct aperture of the phase contrast device for the objective |
| | • Light intensity too low/too high | ➔ Adjust light intensity |
| | • Filter plate is located on the glass panel of the illumination | ➔ Remove filter plate |
| Contrast too strong/edges of objects show distortions | • Aperture diaphragm not set correctly | ➔ Set aperture diaphragm a |
| Microscope image is too bright | • Light intensity too high | ➔ Adjust light intensity |
| | • Aperture diaphragm or field diaphragm were used to control brightness | ➔ Set aperture diaphragm and luminous field diaphragm ➔ Adjust brightness control only via lamp control |
| | • The wrong aperture was set in the phase contrast device | ➔ Select the correct aperture of the phase contrast device for the objective |
| | • Wrong mains adaptor used | ➔ Use the mains adapter supplied |
| Microscope image is too dark | • Light intensity too low | ➔ Adjust light intensity |
| | • Filter plate is located on the glass panel of the illumination | ➔ Remove filter plate |
| | • The specimen between the slide and the cover glass is thicker than necessary | ➔ Microscopy with only as much specimen as necessary |
| | • Dirt on the glass pane of the illumination | ➔ Clean the glass pane carefully. To do this, use compressed air or a fine, lint-free cloth and water/Isopropanol. |
| | • Dirt on or in the condenser | ➔ Carefully clean the lens of the condenser from both sides. To do this, use compressed air or a fine, lint-free cloth and water/Isopropanol. |
| | • Condenser not inserted correctly into the condenser holder | ➔ Push the condenser into the condenser mount as far as possible |
| | • Objective not correctly engaged in the revolving nosepiece | ➔ Turn the revolving nosepiece until it engages correctly |

| FAULT | POSSIBLE CAUSE | REMEDY |
|---|---|--|
| Microscope image is too dark | <ul style="list-style-type: none"> Aperture diaphragm or field diaphragm were used to control brightness | <ul style="list-style-type: none"> → Set aperture diaphragm and luminous field diaphragm → Adjust brightness control only via lamp control |
| | <ul style="list-style-type: none"> The wrong aperture was set in the phase contrast device | <ul style="list-style-type: none"> → Select the correct aperture of the phase contrast device for the objective |
| | <ul style="list-style-type: none"> Wrong mains adaptor used | <ul style="list-style-type: none"> → Use the mains adapter supplied |
| | <ul style="list-style-type: none"> Illumination not switched on or defective | <ul style="list-style-type: none"> → Turn on the illumination → Contact customer service |
| Objects are not visible | <ul style="list-style-type: none"> Objects are too small to be resolved | <ul style="list-style-type: none"> → No error; observe solution |
| | <ul style="list-style-type: none"> The specimen between the slide and the cover glass is thicker than necessary | <ul style="list-style-type: none"> → Microscopy with only as much specimen as necessary |
| Image is orange / green / blue / red | <ul style="list-style-type: none"> Filter on fluorescence unit wrong | <ul style="list-style-type: none"> → Select the correct filter for the fluorescence unit or do not select a filter for bright field measurement |
| Fluorescence unit not working/ Fluorescence microscopy not working | <ul style="list-style-type: none"> Power supply of the fluorescence unit not switched on | <ul style="list-style-type: none"> → Power supply of the fluorescence unit switched on |
| | <ul style="list-style-type: none"> Connection between power supply of fluorescence unit and fluorescence unit is not connected | <ul style="list-style-type: none"> → Connect the connecting cable between the power supply of the fluorescence unit and the fluorescence unit |
| | <ul style="list-style-type: none"> Slider for adjusting the fluorescent LED is neither fully pushed in nor fully pulled out | <ul style="list-style-type: none"> → Either push the slider all the way in or pull it all the way out |
| | <ul style="list-style-type: none"> Filter on fluorescence unit wrong | <ul style="list-style-type: none"> → Select the correct filter of the fluorescence unit |
| | <ul style="list-style-type: none"> Slider for beam splitter of the fluorescence unit inserted | <ul style="list-style-type: none"> → Pull out the slider for the beam splitter of the fluorescence unit completely |
| Incorrect fluorescence wavelength | <ul style="list-style-type: none"> Fluorescence LED adjustment slider not set correctly | <ul style="list-style-type: none"> → Select the correct fluorescent LED using the slider |
| | <ul style="list-style-type: none"> Filter on fluorescence unit wrong | <ul style="list-style-type: none"> → Select the correct filter of the fluorescence unit |
| Image too dark in fluorescence microscopy | <ul style="list-style-type: none"> Aperture diaphragm of the fluorescence unit closed too far | <ul style="list-style-type: none"> → Set aperture diaphragm |
| | <ul style="list-style-type: none"> Light intensity of the fluorescent LED too low | <ul style="list-style-type: none"> → Adjust light intensity |
| | <ul style="list-style-type: none"> Wrong mains adaptor used | <ul style="list-style-type: none"> → Use the mains adapter supplied |
| | <ul style="list-style-type: none"> The specimen between the slide and the cover glass is thicker than necessary | <ul style="list-style-type: none"> → Microscopy with only as much specimen as necessary |
| | <ul style="list-style-type: none"> Objective not correctly engaged in the revolving nosepiece | <ul style="list-style-type: none"> → Turn the revolving nosepiece until it engages correctly |

| FAULT | POSSIBLE CAUSE | REMEDY |
|--|--|--|
| Image too dark in fluorescence microscopy | <ul style="list-style-type: none"> ❖ Slider for adjusting the fluorescent LED is neither fully pushed in nor fully pulled out | <ul style="list-style-type: none"> ➔ Either push the slider all the way in or pull it all the way out |
| | <ul style="list-style-type: none"> ❖ Fluorescence LED switched off or defective | <ul style="list-style-type: none"> ➔ Turn on fluorescent LED ➔ Contact customer service |
| Image too bright in fluorescence microscopy | <ul style="list-style-type: none"> ❖ Aperture diaphragm of the fluorescence unit too wide open | <ul style="list-style-type: none"> ➔ Close the aperture diaphragm further |
| | <ul style="list-style-type: none"> ❖ Light intensity of the fluorescent LED too high | <ul style="list-style-type: none"> ➔ Adjust light intensity |
| | <ul style="list-style-type: none"> ❖ Wrong mains adaptor used | <ul style="list-style-type: none"> ➔ Use the mains adaptor supplied |
| Microscopy image not visible in camera | <ul style="list-style-type: none"> ❖ Slider for switching between binocular and trinocular not set to trinocular | <ul style="list-style-type: none"> ➔ Pull out the slider completely to switch between binocular and trinocular |
| | <ul style="list-style-type: none"> ❖ Camera not switched on/connected | <ul style="list-style-type: none"> ➔ Turn on/connect the camera |
| | <ul style="list-style-type: none"> ❖ All causes, as in No microscope image visible | <ul style="list-style-type: none"> ➔ Take remedial action accordingly <i>No microscope image visible</i> |

➤ **Important note:** Please observe the safety information provided; improper troubleshooting will always invalidate the warranty.

① If the specified remedial measures do not fix the fault, please contact A.KRÜSS Optronic GmbH customer service. You can reach us via the Contact form on the website.