

DENSITY SAMPLES



| Sorted by value | | Alphabetical order | | |
|---------------------------------|-------------------------------|--------------------|--------------------------|---------------------------------|
| ρ [kg/dm ³] | SAMPLES | T [°C] | SAMPLES | ρ [kg/dm ³] |
| 0.621 | Pentane | 20 | Acetic acid | 1.049 |
| 0.640 | Fumaric acid | 20 | Acetone | 0.799 |
| 0.654 | Hexane | 20 | Acetonitrile | 0.783 |
| 0.684 | Heptane | 20 | Acetonylacetone | 0.971 |
| 0.699 | 1-Heptene | 20 | Acetophenone | 1.026 |
| 0.703 | Octane | 20 | Acetyl chloride | 1.103 |
| 0.714 | Diethyl ether | 20 | Acetylacetone | 0.970 |
| 0.718 | 1-Octene | 20 | Acrolein | 0.841 |
| 0.733 | 1-Nonene | 20 | 1-Bromopentane | 1.223 |
| 0.743 | 1-Decene | 20 | Aniline | 1.022 |
| 0.747 | Di-n-propyl ether | 20 | Benzaldehyde | 1.046 |
| 0.752 | 1-Undecene | 20 | Benzene | 0.878 |
| 0.764 | Methylcyclohexane | 20 | Benzyl alcohol | 1.045 |
| 0.779 | Cyclohexane | 20 | Benzyl chloride | 1.098 |
| 0.780 | 1-Pentadecene | 20 | Benzylacetone | 0.989 |
| 0.783 | Acetonitrile | 20 | Bromoform | 2.890 |
| 0.789 | tert-Butanol | 20 | Butyric acid | 0.959 |
| 0.789 | Ethanol | 20 | Caproic acid | 0.929 |
| 0.792 | Methanol | 20 | Caprylic acid | 0.910 |
| 0.799 | Acetone | 20 | Carbon disulfide | 1.174 |
| 0.800 | Methyl isobutyl ketone (MIBK) | 20 | Carvacrol | 0.976 |
| 0.802 | Isobutanol | 20 | Caustic hypochlorite | 1.220 |
| 0.804 | 1-Propanol | 20 | Chlorobenzene | 1.107 |
| 0.805 | Cyclopentadiene | 20 | Chlorocyclohexane | 1.000 |
| 0.805 | Methyl ethyl ketone (MEK) | 20 | Chloroform | 1.489 |
| 0.810 | 1-Butanol | 20 | cis-1,2-Dichloroethylene | 1.282 |
| 0.810 | Kerosene | 20 | cis-1,2-Dibromoethylene | 2.246 |
| 0.811 | Cyclohexene | 20 | Citral | 0.893 |

DENSITY SAMPLES



| Sorted by value | | | Alphabetical order | |
|---------------------------------|---------------------------------------|-----------|-----------------------------------|---------------------------------|
| ρ [kg/dm ³] | SAMPLES | T [°C] | SAMPLES | ρ [kg/dm ³] |
| 0.814 | 3-Heptanone | 20 | Crotonaldehyde | 0.856 |
| 0.816 | 1-Pentanol | 20 | Cyclohexane | 0.779 |
| 0.820 | 1-Hexanol | 20 | Cyclohexanol | 0.962 |
| 0.823 | 1-Heptanol | 20 | Cyclohexanone | 0.949 |
| 0.826 | Dimethylamine solution 60 % (DMA) | 20 | Cyclohexene | 0.811 |
| 0.827 | 1-Octanol | 20 | Cyclohexylamine | 0.896 |
| 0.828 | 1-Nonanol | 20 | Cyclopentadiene | 0.805 |
| 0.829 | 1-Decanol | 20 | Dichloromethane | 1.336 |
| 0.841 | Acrolein | 20 | Diethyl ether | 0.714 |
| 0.850 | Mesityl oxide | 20 | Diethylaniline | 0.934 |
| 0.853 | Lemon oil | 20 | Dimethylacetamide 90 % | 0.940 |
| 0.856 | Crotonaldehyde | 20 | Dimethylamine solution 60 % (DMA) | 0.826 |
| 0.860 | Piperidine | 20 | Dimethylaniline | 0.956 |
| 0.863 | Linalool | 20 | Dimethylformamide (DMF) | 0.948 |
| 0.863 | Mesitylene | 20 | Di-n-propyl ether | 0.747 |
| 0.863 | p-Xylene | 20 | 1-Decene | 0.743 |
| 0.866 | 1-Chlorodecane | 20 | Enanthic acid | 0.922 |
| 0.868 | Ethylbenzene | 20 | Ethanol | 0.789 |
| 0.868 | 1,2,4-Trimethylbenzene (Pseudocumene) | 20 | Ethyl acetate | 0.900 |
| 0.870 | Toluene | 20 | Ethyl iodid | 1.940 |
| 0.871 | m-Xylene | 20 | Ethylbenzene | 0.868 |
| 0.872 | 1-Chlorohexane | 20 | Ethylen bromide | 2.180 |
| 0.872 | 1-Chlorooctane | 20 | Ethylenglycole | 1.115 |
| 0.878 | Benzene | 20 | Formanide | 1.139 |
| 0.878 | Isopropylbenzene (Cumene) | 20 | Formic acid | 1.212 |
| 0.884 | 1-Chlorobutane | 20 | Fumaric acid | 0.640 |
| 0.887 | Hemellitol | 20 | Glycerol | 1.261 |
| 0.891 | Propyl acetate | 20 | Hemellitol | 0.887 |



DENSITY SAMPLES



| Sorted by value | | | Alphabetical order | |
|---------------------------------|---------------------------------|-----------|---------------------------------|---------------------------------|
| ρ [kg/dm ³] | SAMPLES | T [°C] | SAMPLES | ρ [kg/dm ³] |
| 0.893 | Citral | 20 | Heptane | 0.684 |
| 0.896 | Cyclohexylamine | 20 | Hexamethylenediamine adipinat | 1.201 |
| 0.900 | Ethyl acetate | 20 | Hexane | 0.654 |
| 0.900 | Methylamine solution 40 % (MMA) | 20 | Indene | 0.998 |
| 0.910 | Caprylic acid | 20 | Iodobenzene | 1.830 |
| 0.922 | Enanthic acid | 20 | Iron(II) sulfate heptahydrate | 1.900 |
| 0.929 | Caproic acid | 20 | Isobutanol | 0.802 |
| 0.934 | Diethylaniline | 20 | Isopropylbenzene (Cumene) | 0.878 |
| 0.940 | Dimethylacetamide 90 % | 20 | Kerosene | 0.810 |
| 0.942 | Vinyl acetate (VyAc) | 20 | Lemon oil | 0.853 |
| 0.948 | Dimethylformamide (DMF) | 20 | Linalool | 0.863 |
| 0.949 | Cyclohexanone | 20 | Lithium chloride | 2.068 |
| 0.956 | Dimethylaniline | 20 | Mercury | 13.595 |
| 0.959 | Butyric acid | 20 | Mesityl oxide | 0.850 |
| 0.962 | Cyclohexanol | 20 | Mesitylene | 0.863 |
| 0.967 | Tetralin | 20 | Methanol | 0.792 |
| 0.970 | Acetylacetone | 20 | Methyl diethanolamin (MDEA) | 1.040 |
| 0.971 | Acetylacetone | 20 | Methyl ethyl ketone (MEK) | 0.805 |
| 0.976 | Carvacrol | 20 | Methyl iodide | 2.279 |
| 0.984 | N-Methylaniline | 20 | Methyl isobutyl ketone (MIBK) | 0.800 |
| 0.989 | Benzylacetone | 20 | Methylamine solution 40 % (MMA) | 0.900 |
| 0.992 | Propionic acid | 20 | Methylcyclohexane | 0.764 |
| 0.994 | Paraldehyde | 20 | m-Toluidine | 0.998 |
| 0.998 | Indene | 20 | m-Xylene | 0.871 |
| 0.998 | m-Toluidine | 20 | N-Benzyl-N-ethylaniline | 1.029 |
| 0.998 | o-Toluidine | 20 | Nicotine | 1.009 |
| 0.998 | Water | 20 | Nitro benzene | 1.207 |
| 1.000 | Chlorocyclohexane | 20 | Nitro methane | 1.139 |



DENSITY SAMPLES



| Sorted by value | | T [°C] | Alphabetical order | |
|---------------------------------|-------------------------------------|-----------|--|---------------------------------|
| ρ [kg/dm ³] | SAMPLES | | SAMPLES | ρ [kg/dm ³] |
| 1.009 | Propiophenone (Ethyl phenyl ketone) | 20 | 3-Nitrotoluene | 1.157 |
| 1.009 | Nicotine | 20 | 2-Nitrotoluene | 1.163 |
| 1.022 | Aniline | 20 | N-Methylaniline | 0.984 |
| 1.026 | Acetophenone | 20 | Octane | 0.703 |
| 1.029 | N-Benzyl-N-ethylaniline | 20 | 1-Bromooctane | 1.166 |
| 1.038 | 1,4-Dioxane | 20 | o-Toluidine | 0.998 |
| 1.040 | Methyl diethanolamin (MDEA) | 20 | Paraldehyde | 0.994 |
| 1.045 | Benzyl alcohol | 20 | Pentachloroethane | 1.672 |
| 1.046 | Benzaldehyde | 20 | Pentane | 0.621 |
| 1.049 | Acetic acid | 20 | Phenylhydrazine | 1.098 |
| 1.064 | Thiophene | 20 | Phthalic anhydride | 1.527 |
| 1.066 | 4-Chlorotoluene | 20 | Piperidine | 0.860 |
| 1.069 | Quinaldine | 20 | Propionic acid | 0.992 |
| 1.070 | 3-Chlorotoluene | 20 | Propiophenone (Ethyl phenyl ketone) | 1.009 |
| 1.093 | Quinoline | 20 | Propyl acetate | 0.891 |
| 1.098 | Benzyl chloride | 20 | p-Xylene | 0.863 |
| 1.098 | Phenylhydrazine | 20 | Pyruvic acid | 1.267 |
| 1.103 | Acetyl chloride | 20 | Quinaldine | 1.069 |
| 1.107 | Chlorbenzene | 20 | Quinoline | 1.093 |
| 1.115 | Ethylenglycole | 20 | Sodium hydroxide solution (Caustic soda) | 1.430 |
| 1.139 | Formanide | 20 | Sulfuric acid 90 % | 1.810 |
| 1.139 | Nitro methane | 20 | tert-Butanol | 0.789 |
| 1.157 | 3-Nitrotoluene | 20 | Tetrabromomethane | 2.963 |
| 1.163 | 2-Nitrotoluene | 20 | Tetrachloroethylene | 1.614 |
| 1.166 | 1-Bromooctane | 20 | Tetralin | 0.967 |
| 1.174 | Carbon disulfide | 20 | Thiophene | 1.064 |
| 1.201 | Hexamethylenediamine adipinat | 20 | Toluene | 0.870 |
| 1.207 | Nitro benzene | 20 | trans-1,2-Dibromoethylene | 2.231 |



DENSITY SAMPLES



| Sorted by value | | Alphabetical order | | |
|---------------------------------|--|--------------------|---------------------------------------|---------------------------------|
| ρ [kg/dm ³] | SAMPLES | T [°C] | SAMPLES | ρ [kg/dm ³] |
| 1.212 | Formic acid | 20 | trans-1,2-Dichloroethylene | 1.257 |
| 1.220 | Caustic hypochlorite | 20 | Tribromoacetaldehyde (Bromal) | 2.550 |
| 1.223 | 1-Bromopentane | 20 | Vinyl acetate (VyAc) | 0.942 |
| 1.253 | 1,2-Dichlorethane | 20 | Water | 0.998 |
| 1.257 | trans-1,2-Dichloroethylene | 20 | 1,1,2,2-Tetrabromoethane | 2.970 |
| 1.261 | Glycerol | 20 | 1,2,4-Trichlorobenzene | 1.477 |
| 1.267 | Pyruvic acid | 20 | 1,2,4-Trimethylbenzene (Pseudocumene) | 0.868 |
| 1.275 | 1-Bromobutane | 20 | 1,2-Dichlorethane | 1.253 |
| 1.282 | cis-1,2-Dichloroethylene | 20 | 1,2-Dichlorobenzene | 1.305 |
| 1.296 | 2-Nitroethanol | 20 | 1,4-Dioxane | 1.038 |
| 1.305 | 1,2-Dichlorobenzene | 20 | 1-Bromobutane | 1.275 |
| 1.336 | Dichloromethane | 20 | 1-Bromonaphthalene | 1.487 |
| 1.430 | Sodium hydroxide solution (Caustic soda) | 20 | 1-Butanol | 0.810 |
| 1.441 | 1-Iodoheptane | 20 | 1-Chlorobutane | 0.884 |
| 1.477 | 1,2,4-Trichlorobenzene | 20 | 1-Chlorodecane | 0.866 |
| 1.487 | 1-Bromonaphthalene | 20 | 1-Chlorohexane | 0.872 |
| 1.489 | Chloroform | 20 | 1-Chlorooctane | 0.872 |
| 1.527 | Phthalic anhydride | 20 | 1-Decanol | 0.829 |
| 1.614 | 1-Iodobutane | 20 | 1-Heptanol | 0.823 |
| 1.614 | Tetrachloroethylene | 20 | 1-Heptene | 0.699 |
| 1.672 | Pentachloroethane | 20 | 1-Hexanol | 0.820 |
| 1.747 | 1-Iodopropane | 20 | 1-Iodobutane | 1.614 |
| 1.810 | Sulfuric acid 90 % | 20 | 1-Iodoheptane | 1.441 |
| 1.830 | Iodobenzene | 20 | 1-Iodopropane | 1.747 |
| 1.900 | Iron(II) sulfate heptahydrate | 20 | 1-Nonanol | 0.828 |
| 1.940 | Ethyl iodid | 20 | 1-Nonene | 0.733 |
| 2.068 | Lithium chloride | 20 | 1-Octanol | 0.827 |
| 2.180 | Ethylen bromide | 20 | 1-Octene | 0.718 |



DENSITY SAMPLES



| Sorted by value | | | Alphabetical order | |
|---------------------------------|-------------------------------|-----------|--------------------|---------------------------------|
| ρ [kg/dm ³] | SAMPLES | T [°C] | SAMPLES | ρ [kg/dm ³] |
| 2.231 | trans-1,2-Dibromoethylene | 20 | 1-Pentadecene | 0.780 |
| 2.246 | cis-1,2-Dibromoethylene | 20 | 1-Pentanol | 0.816 |
| 2.279 | Methyl iodide | 20 | 1-Propanol | 0.804 |
| 2.550 | Tribromoacetaldehyde (Bromal) | 20 | 1-Undecene | 0.752 |
| 2.890 | Bromoform | 20 | 2-Nitroethanol | 1.296 |
| 2.970 | 1,1,2,2-Tetrabromoethane | 20 | 3-Chlorotoluene | 1.070 |
| 2.963 | Tetrabromomethane | 20 | 3-Heptanone | 0.814 |
| 13.595 | Mercury | 20 | 4-Chlorotoluene | 1.066 |

SAMPLES

The table presents an overview of samples that our customers typically measure with our density meters.
The values should not be interpreted as specifications.

A.KRÜSS Optronic GmbH
Alsterdorfer Straße 276–278
22297 Hamburg | Germany

Phone +49 40 514317-0
Fax +49 40 514317-60

E-Mail info@kruess.com
Web www.kruess.com

